




tANCHOR

Vector kit



For presenting peptides on the surface of mammalian cells

ATG:biosynthetics catalog no. TA-CDA-ASA1-1

User Manual Version V2.3

2018-08-01

100111ATG: 
ATG:biosynthetics
Your Partner in Synthetic Biology

ATG:biosynthetics GmbH
Weberstraße 40 • 79249 Merzhausen • Germany
ms@atg-biosynthetics.com • www.atg-biosynthetics.com

This kit is for research only!

Not for diagnostic purposes!

This product requires knowledge of basic molecular biological techniques

Notice for purchased products

Products from ATG:biosynthetics GmbH are to be used for research purposes only. Products are not made for use *in vitro* diagnostic applications, therapeutics, or in humans. ATG:biosynthetics does not allow to transfer rights to third parties or reselling products without permission or selling modified products. Further, products from ATG:biosynthetics can't be used to provide services to third parties without the permission of ATG:biosynthetics, GmbH, Merzhausen, Germany on a license basis.

Limited Use Label License No: TA-1 ATG:biosynthetics GmbH

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) inclusion of the product as part of another product; (3) use of the product or its components to provide a service, information, or data; (4) use of the product or its components for *ex vivo* or *in vivo* therapeutic, diagnostic, prophylactic or other unauthorized commercial purposes; (5) use in foods, drugs, devices or cosmetics of any kind, or for consumption by or use in connection with or administration or application to humans or animals; (6) resale of the product or its components, whether or not such product or its components are resold for use in research, or (7) any use other than research use. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. ATG will not assert a claim against the buyer of infringement of patents owned or controlled by ATG which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser does not agree with the terms or limitations of this limited use agreement, the customer can return the product at their own costs but for a full refund (excluding shipping, tax and handling fees). For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, ATG:biosynthetics GmbH, D-79249 Merzhausen, Weberstr. 40, Germany or



outlicensing@atg-biosynthetics.com. However, in accepting ATG's material transfer agreement, buyers are not discouraged from disseminating the material to their colleagues if the recipient is willing to sign this agreement prior to transfer of the material combined with its regular purchase from ATG. *tANCHOR* - vectors are for research purposes only. By signing this order form and on receipt of the *tANCHOR* vectors you agree not to distribute *tANCHOR* vectors yourself and that you will only use the vectors in your individual lab environment. For all other purposes a specific and separate license needs to be purchased. Third parties who come into possession of the *tANCHOR* vectors through your distributing them need to sign the same agreement, need to formally order the vectors at ATG and reimburse ATG for permission to use them.

Data Safety and Confidentiality Statement:

ATG:biosynthetics declares that we will keep your intellectual property about your genes/ sequences in strict confidentiality and we will NOT, in any manner, disclose, develop and utilize your genes/ sequences. We will destroy all the related testing data and results after at least one month of storing. On request we will provide a secure web space for mutual data exchange.

Security Statement:

Please note that ATG like other gene synthesis companies in compliance with bio-security relevant self-control has subjected itself to a "code of conduct" which is emerging to be an international quality management standard of the world wide community of Synthetic Biologists. Sequences are checked for their potential of being a bio-hazard or their potential for misuse in bio-terrorism. If you know that your gene has any such potential, please cooperate with the manufacturer or distributor, e.g. ATG.

By buying this product you are accepting the license conditions and agree.

Otherwise you will send immediately back all items of the *tANCHOR* kit without keeping items back or share it with third parties.



Table of contents

1. Product information	3
2. Introduction	4
2.1 Application	4
3. Methods	6
3.1. Using the right <i>tANCHOR</i> vector.....	6
3.2. Cloning of epitope sequences in the <i>tANCHOR</i> system	6
3.3. Vector features	7
3.4. Vector maps.....	7
3.5. Cloning diagram	9
3.6. Transfection of human cells	9
3.7. Verification of protein expression	10
3.8. Detection of displayed epitopes	10
4. Troubleshooting	11
5. References	12



1. Product information

Content: 5 µg of each *t*ANCHOR vector (lyophilized in TE buffer)

Storage: Long term: lyophilized vectors should be stored at -20°C
Short term: TE buffer 4°C or in distilled water at -20°C

Product characterization: *t*ANCHOR vectors can be characterized by restriction enzymes showing the correct following banding pattern

Vector	Enzyme	Fragment length in bp
<i>t</i> ANCHOR-CD63	<i>NheI/ClaI</i>	621 and 3535
<i>t</i> ANCHOR-CD63	<i>XhoI/EcoRI</i>	1098 and 3058
<i>t</i> ANCHOR-CD82	<i>NheI/ClaI</i>	666 and 3535
<i>t</i> ANCHOR-CD82	<i>XhoI/EcoRI</i>	1113 and 3088

License number: Each purchased *t*ANCHOR kit will be shipped with a unique license number. The license number is stamped on the *t*ANCHOR kit package. Please keep the license number (customer#) for your records, it will enable you to receive quick and personalized support. Please consult us if you encounter any problems during your applications with the *t*ANCHOR vectors. We will help you as quickly as possible.

2. Introduction

The *tANCHOR* technology provides an innovative and versatile tool for displaying epitopes on the surface of e.g. human cells which was developed by Dr. Daniel Ivanusic. This system is based on the use of tetraspanin anchors (*tANCHOR*). Transmembrane anchors derived from a bundle of 4 tetraspanin transmembrane domains assembled as a type III integral membrane protein are highly efficient to transport epitopes to the surface of human cells. The main characteristics of tetraspanin protein family consist of the presence of a small and a large extracellular loop (SEL and LEL) which are flanked by four transmembrane domains. The epitope sequence for displaying on the cell surface is placed instead of the large extracellular loop of the tetraspanin full-length sequence and is flanked by the restriction sites *EcoRI* and *EcoRV* (Fig. 1).




Fig. 1

2.1 Application

During its development the *tANCHOR* system was verified and tested for functionality by displaying epitopes 2F5-4E10 derived from HIV-1, V5-6xHis, relevant epitopes for checkpoint targeted therapies of human PD-1 and CD20 on the cell surface. Confocal analysis of *tANCHOR* proteins fused to mCherry (Fig. 2) showed correctly expression on the surface of human cells (HEK293T, HeLa). The extracellular orientation was tested with anti-V5 in an ELISA based assay and showed that epitopes are highly extracellular accessible. Expressed *tANCHOR* proteins can be used for bioanalytical methods (e.g. ELISA), protein-

protein-interaction studies on the cell surface or as antigen expression system for immunization studies.





Fig. 2

We recommend that the constructs should be fused C-terminal with a mCherry fusion-tag. This fusion will significantly reduce the time needed for localization verification of expressed proteins. Delivered *tANCHOR* vectors do not contain fusion of a mCherry sequence.

tANCHOR technology provides an innovative tool for displaying epitopes on the surface of human cells. This enables you to improve your ELISA methods to search for antibodies binding studies under native conditions on your target (**Fig. 3**).



1. Seed cells
2. Transfect and incubate
3. Antibody incubation
4. Washing steps
5. HRP detection (TMB)

Fig. 3

The benefits are:

- Maximize your search for antibodies under native conditions
- No need to coat plates with expensive peptides
- High throughput screening: seed and transfect cells – incubate your antibodies
- Fast and reliable results

3. Methods

3.1. Using the right *tANCHOR* vector

The *tANCHOR* FULLkit contains two *tANCHOR* ASSAY-vectors *tANCHOR*-CD63 and *tANCHOR*-CD82 and 6 control vectors. We recommended start with *ptANCHOR*-CD82 if you are using HeLa as target cells. If you are using cells where endogenous CD63 is highly expressed, then use *tANCHOR*-CD63.

3.2. Cloning of epitope sequences in the *tANCHOR* system

Upon receiving the vial resuspend the lyophilized *tANCHOR* vectors in 50 μ l sterile water to make a 0.1 μ g/ μ l stock solution. Vortex for 1 minute and store the stock solution at -20°C . Use this stock solution to transform an *E. coli* strain such as TOP10 or DH5alpha. If you are using other *E. coli* strains, ensure that they are recombination (*recA*) and endonuclease A deficient (*endA*). The *tANCHOR* vector system is a fusion protein expression vector requiring that you have to clone your epitope sequence in frame between restriction sites *EcoRI* and *EcoRV* with the N- and the C-terminal part of the tetraspanin transmembrane sequence (Fig. 4). See sequence map and cloning diagram for designing your cloning strategy (Fig. 5-7). Please consult us at support@atg-biosynthetics.com if you have further questions. The sequences of *tANCHOR* vectors are delivered with the *tANCHOR* kit or are available upon request.




Fig. 4

3.3. Vector features

The *tANCHOR* vectors contains the following feature elements.

Feature	Restriction sites	Function	Position	
			ptANCHOR-CD63	ptANCHOR-CD82
<i>NcoI-X-NcoI</i>	CC AT GG	FLAG-tag	732-740	732-740
<i>BstXI</i>	CCAN NNNN NTGG	Alternative site for introduction of 5'-tags	762-773	762-773
<i>NcoI-EcoRI</i>	CC AT GG – G AATT C	TM1-3	761-1105	762-1115
<i>EcoRI-EcoRV</i>	G AATT C – GAT ATC	V5-6xHis-tag control epitopes	1100-1201	1115-1216
<i>EcoRV-ClaI</i>	GAT ATC – AT CG AT	TM4	1196-1339	1211-1384
<i>ClaI-PmeI</i>	AT CG AT – GTT AAAC	c-myc-Tag control epitope	1334-1377	1384-1422
<i>PmeI</i>	GTT AAAC	Remove empty vectors	1370-1377	1415-1422

Table 1

3.4. Vector maps




Fig. 5




Fig. 6

In order to introduce your epitope sequence in the *tANCHOR* vector create a fragment by PCR or gene synthesis with the flanked by the sequence of appropriate restriction sites GAATTC (*EcoRI*) and GATATC (*EcoRV*). Clone these fragments using standard cloning procedures. Note: in order to cut effectively 4-6 nucleotides should be added on either side of the recognition site. Detailed information is available from the restriction enzyme supplier like on the homepage of New England Biolabs (Cleavage Close to the End of DNA Fragments). After incubation of your ligation batch over night at 16°C transform into chemically competent *E. coli* cells, plate cells on LB agar plates containing 100 µg/ml Ampicillin, pick clones and check recovered plasmids for the presence of the insert of interest using restriction enzymes and confirm sequence by Sanger DNA sequencing. In order to keep the rate of empty vectors low you can incubate your ligation batch with *PmeI*.


Use the following sequence primers to confirm your construct:

	Primer sequence 5'-3'
Seq1_for	GGATAGCGGTTTGA CT CACG
Seq1_rev	GCACAGTCGAGGCTTTATGT

3.5. Cloning diagram

Carefully inspect your cloning strategy with the following map. Your epitope sequence must be *in frame* with the *tANCHOR* sequence as shown in the cloning diagram. The control vector contains a V5-6xHis insert which is in frame between the *EcoRI* and *EcoRV* restriction sites and serves as an example for the cloning strategy. We recommend inserting your fragment between the *EcoRI/EcoRV* restriction sites. The V5 or 6xHis epitope can be fused with your epitope sequence by inserting your sequences between *EcoRI/AgeI* or *AgeI/EcoRV*.

tANCHOR-CD63



tANCHOR-CD82




Fig. 7

3.6. Transfection of human cells

This section is provided for your convenience to increase research success. More details are provided in the references.

Once the correct sequence of your clone is confirmed you will have to prepare plasmid DNA for transfection of cells. Isolated plasmid DNA for transfection into human cells have to be highly pure and containing low endotoxin (LPS) contaminations. LPS contaminations are able to kill cells or decrease



transfection efficiency. Impurities with salt will in addition lead to decrease your transfection efficiency. We recommended to prepare plasmid DNA using the EndoFree MaxiPrep kit from Qiagen, Germany. For transfection of HEK293T or HeLa cells we recommended MetafectenePro from Biontech, Germany for high transfection efficiency of *tANCHOR* vectors according to manufacturer's information. Jurkat cells (immortalized cell line of T lymphocyte cells) were also tested using AMAXA transfection kit with the nucleofection program S-10. Insertion of mCherry (see vector map for cloning strategy) instead of the cMyc tag will easy display ease display transfection efficiency. The use of a transfection control e.g. pmCherry-N1 from Clontech, Germany is useful.

3.7. Verification of protein expression

Verify your protein expression by Western blot analysis from total cell extracts using either an antibody directed against your inserted epitope of interest or an antibody directed against the N-terminal FLAG-tag or C-terminal cMyc-tag of the fusion protein. After confirming a positive protein expression proceed with the protein localization analysis.

3.8. Detection of displayed epitopes

Basic *tANCHOR* vectors does not include a fluorescence reporter protein e.g. mCherry therefore displayed epitopes have to be visualized by indirect immunofluorescence technique. We recommended to detect the V5 epitope in the basic vector with an anti-V5 antibody and a secondary fluorescently conjugated antibody. Instead of detecting the extracellular orientation by indirect immunofluorescence there is also the possibility to test by standard ELISA-based methods. For this application seed cells in a 96 well plate, transfect and test the indicated presented peptide by an antibody recognizing an epitope within the peptide (**Fig. 3**). The bound antibody can be detected an HRP-conjugated secondary antibody. Use a control antibody e.g. anti-cMyc for excluding intracellular or unspecific antibody binding.

4. Troubleshooting

Issue	Possible reason	Solution
No protein expression (western blot analysis)	Low transfection efficiency	Improve transfection efficiency
	Low DNA quality	Prepare new plasmid DNA batch
	Epitope sequence is not in frame	Check sequence
No peptide presentation on the cell surface (method: indirect immunofluorescence)	Failure during performing indirect immunofluorescence techniques	Use other primary or fluorescently conjugated secondary antibody
		If issue still exist clone mCherry sequence C-terminally
No peptide presentation on the cell surface (method: direct immunofluorescence)	Failure during transfection	Improve transfection efficiency
	mCherry is not in frame	Check sequence
Control antibody binds on the cells (method: ELISA based)	Antibody concentration too high	Make dilutions of your control antibody
	Cells dried during washing steps → plasma membrane breaks and intracellular epitopes becoming accessible for control antibodies	Ensure that cells never dried while performing ELISA steps



5. References

Background

Boucheix, C., and Rubinstein, E. (2001). Tetraspanins. *Cell Mol Life Sci* 58, 1189-1205.

Kovalenko, O.V., Metcalf, D.G., DeGrado, W.F., and Hemler, M.E. (2005). Structural organization and interactions of transmembrane domains in tetraspanin proteins. *BMC Struct Biol* 5, 11.

Stipp, C.S., Kolesnikova, T.V., and Hemler, M.E. (2003). Functional domains in tetraspanin proteins. *Trends Biochem Sci* 28, 106-112.

Ivanusic, D (2016). System for presenting peptides on the cell surface. European Patent No. EP20160710935, Munich, Germany: European Patent Office.

Methods

Green, M.R. and Sambrook, J. (2012). *Molecular Cloning: A Laboratory Manual*, 4th Edition New York: Cold Spring Harbor Laboratory Press.

Matsumura I. Why Johnny can't clone: Common pitfalls and not so common solutions. *BioTechniques*. 2015;59(3):IV-XIII.

Ivanusic D, Pietsch H, Konig J, Denner J. Absence of IL-10 production by human PBMCs co-cultivated with human cells expressing or secreting retroviral immunosuppressive domains. *PLoS One*. 2018;13(7):e0200570.

Ivanusic D, Eschricht M, Denner J. Investigation of membrane protein-protein interactions using correlative FRET-PLA *BioTechniques*. 2014;57:188-98.

Donaldson JG. Immunofluorescence Staining. *Curr Protoc Cell Biol*. 2015;69:4 3 1-7.

Chicaybam L, Sodre AL, Curzio BA, Bonamino MH. An Efficient Low Cost Method for Gene Transfer to T Lymphocytes. *PLoS ONE*. 2013;8(3):e60298.

Ivanusic D, Madela K, Denner J. Easy and cost-effective stable positioning of suspension cells during live-cell imaging. 2017. 2017.

Mahmood T, Yang P-C. Western Blot: Technique, Theory, and Trouble Shooting. *North American Journal of Medical Sciences*. 2012;4(9):429-34.